





SIDE VIEW


FRONT VIEW



Rear Underbrace

Approval _____

Date _____

	Approvals Drawn By J. Lamb Approved By		Date 4/20/11		GAT-TBL-SNK-304036-1-BS-SS-ST
	TITLE 30"Dx40"Wx36"H SS Lab Sink Tbl. 15"Wx15"Dx10"H Sink Faucet & Drain			Material: SS 304 #3	
<small>PROPRIETARY THIS DOCUMENT AND THE DATA DISCLOSED HEREIN OR HERewith IS NOT TO BE REPRODUCED, USED OR OTHERWISE DISCLOSED IN WHOLE OR IN PART TO ANYONE WITHOUT THE WRITTEN AUTHORIZATION FROM G2 AUTOMATED TECHNOLOGIES, LLC</small>					
SIZE N/A		CAGE CODE		DWG NO GAT-TBL-SNK-2028-OVERALL	
SCALE		REV		SHEET	

A4

